

Lecture Bottle Stand



The 518-1616 Lecture Bottle Stand is designed to safely secure a lecture bottle or other 2" diameter gas cylinder to the workbench. These PVC constructed units are rugged and chemically resistant. The base is predrilled with four holes for permanent installation.

Materials and Specifications

- *Material*
PVC
- *Size*
8" x 8" x 5" (20.3 cm x 20.3 cm x 12.7 cm)
- *Weight*
1.02 lbs. (.46 kg)

Ordering Information

Stock Number	Description
518-1616	Lecture Bottle Stand

Lecture Bottle Wall Mount



The 518-1617 wall mount bracket is designed to secure a lecture bottle for wall or bench mounting. These PVC constructed units are rugged and chemically resistant. Predrilled holes in the base and rear panel facilitate mounting.

Materials and Specifications

- *Material*
PVC
- *Size (518-1617)*
8" x 8" x 5" (20.3 cm x 20.3 cm x 12.7 cm)
- *Size (518-1618)*
11" x 4" x 7" (27.9 cm x 10.2 cm x 17.8 cm)
- *Weight (518-1617)*
1.02 lbs. (.46 kg)
- *Weight (518-1618)*
4.0 lbs. (1.8 kg)

Ordering Information

Stock Number	Description
518-1617	Lecture Bottle Wall Mount (capacity: one lecture bottle)
518-1618	Lecture Bottle Wall Mount (capacity: three lecture bottles)